

LiteSafe® Sample Storage Products

Why Put Your Research at Risk!

Argos offers a complete family of products to keep your samples, reagents and solutions safe. Unlike amber tinted tubes, these products are completely opaque, therefore protecting even the most light sensitive compounds. Argos LiteSafe® tubes eliminate the inconvenience of fitting tubes with aluminum foil to protect samples. Ideal for use with fluorescent probes or ethidium bromide.

- Manufactured from polypropylene; autoclavable to 121°C
- Suitable for freezing to -100°C
- 100-Place freezer box accommodates 0.5 - 2.0mL micro-tubes, as well as most brands of cryogenic vials. Includes alpha-numeric grid.
- Tubes are RNase, DNase free
- Frosted writing surface on tubes for sample identification
- 15 and 50mL tubes are polypropylene with HDPE caps

Ordering Information

Cat. No.	Description	Qty.
T7456-001	0.5mL LiteSafe microcentrifuge tube	500/pack
T7100BK	1.5mL LiteSafe microcentrifuge tube	500/pack
T2081A	5mL LiteSafe microcentrifuge tube	200/pack
TB1500	15mL LiteSafe centrifuge tube	50/pack
TB5000	50mL LiteSafe centrifuge tube	50/pack
R3121	100-place LiteSafe storage box	5/pack



Motorized Pestle Mixer

The Argos Pestle Mixer is a lightweight battery powered instrument that can be used with disposable pestles, mixers and micro tubes. Cell tissue and pellets are completely homogenized within seconds by means of a high speed vortexing action. Greatly reduces hand fatigue. Comes complete with 2 AA batteries and pestle/mixer adapter.

Pestles & Tubes

Designed to homogenize cells and plant tissue or to re-suspend protein and DNA pellets. Argos pestles and tubes are DNase / RNase and pyrogen free and individually wrapped to guarantee purity. With their grooved grip handle tailored to fit with our microtube, it is easy to turn and does not slip even when wearing gloves. Autoclavable. Pestle length: 7.5 cm

Ordering Information

Cat. No.	Description	Qty.
A0001	Pestle mixer	1/each
A0002	Pestle mixer replacement rubber adapter (mixer to pestle connector)	1/each
P7338-901	0.5mL Pestle	100/pack
P7339-901	1.5mL Pestle	100/pack
P7811-904	1.5mL Microtube	100/pack
P8185-904	0.5mL Microtube	100/pack
P9950-901	1.5mL Microtube and pestle (1 pestle & 1 microtube per individual bag)	100/pack
P9951-901	0.5mL Microtube and pestle (1 pestle & 1 microtube per individual bag)	100/pack

