

2026 PROMOTION

1 JANUARY - 31 DECEMBER

LAB EQUIPMENTS & CONSUMABLE



LINE Official Account

AS ONE

HOT SALE

18,500.-

Auto Dry Desiccator
Item no. 1-5486-21, 38L
Model OH-3S (220V)

HOT SALE

25,500.-

Auto Dry Desiccator
Item no. 1-5488-22, 38L
Model UVOH-520SA (220V)

7,900.-

Desiccator
Item no. 1-001-01, 38L
Model LH

BIG Size!!

48,000.-

Auto Dry Desiccator
Item no. 1-5504-81
Model SP-BGFN-P (100V)

Humidity can be set **NEW ARRIVAL**

20,500.-

Auto Dry Desiccator
Item no. 4-6388-01
Model ZL-OH-3S (100V)

Humidity can be set **NEW ARRIVAL**

22,700.-

Auto Dry Desiccator
Item no. 4-6389-01
Model ZL-OH-3K,
Light-shielding PVC (100V)

32,000.-

Cool Incubator
Item no. 2-926-21
Temp: 3 to 45 °C, 28 L (220V)

120,000.-

Vacuum Dryer
Item no. 1-2186-12, Model AVO-250V
Temp: 40-120 °C, 15L (100V)
*มีปั๊มสุญญากาศ Vacuum pump

68,000.-

Cooling Water Circulator
Item no. 1-5469-41, Model LTC-450
Temp: -20 to +120 °C, 4.5L (100V)

79,900.-

Constant Temperature Stirrer Bath
Item no. 1-3903-2, Model KS-1X
Temp: R.T. + 5 to 80 °C (100V)
Rotation speed: 100 to 1500 rpm
Stirring capacity: 300 mL Beaker x 8 pc

FREE

31,500.-

Incubation Cover
Item no. 1-7605-01
Temp: Room temperature + 5 - 45°C
Inner size 422 x 372 x 346mm (100V)

29,900.-

Clean Booth AIG Type
Item no. 3-4055-81
dust collection efficiency: 99.97% or more (0.3 μm particles) (220V)

HOT SALE

26,900.-

Magnetic Stirrer 4 Position
Item no. 1-4604-42
Speed (rpm): MAX.1500
Stirring capacity (mL): 50 ~ 2000 x 4 (220V)

HOT SALE

38,500.-

Magnetic Stirrer 6 Position
Item no. 1-4605-42
Speed (rpm): MAX.1500
Stirring capacity (mL): 50 ~ 2000 x 6 (220V)

BIG Size!!

48,000.-

Hot Plate
Item no. 2-644-81, Model HPR-4030G
Highest temperature (°C): 300
Plate size (mm): 400 x 300
Plate material: Aluminum + ceramic coating (220V)

36,000.-

SPE Vacuum Manifold
Item no. 4-2854-01
number of ports : 12
material: Vacuum Chamber/Glass, Lid/Marble, Valve/PP

39,900.-

Moisture activity
Item no. 4-2453-03, Model MD-AW3
measurement range: 0 to 1.00 Aw
measuring method: Electrostatic capacity Humidity Sensor
999 records CSV files can be retrieved using a USB cable.

HOT SALE

61,000.-

UV Lamp Viewing Cabinet with UV Lamp
Item no. 1-5479-22+1-5479-09
Standard electrophoresis gels and TLC plates (200~200 mm) can be observed.
365nm, 254nm switch set (100V)

HOT SALE

23,900.-

Gel Electrophoretic
Item no. 1-5484-01
Output: 18 to 135 V
Timer: 0-99 min (220V)

HOT SALE

2,900.-

Waterproof pH Meter
Item no. 2-1287-22
Measurement range (pH): 0.0~14.0
Waterproof (IP57 compliant)

18,500.-

Waterproof Portable Conductivity Meter
Item no. 1-2814-01
Conductivity, TDS, salinity, resistivity, Temp
Measurement range: Conductivity 0.00 μS/cm - 199.9mS/cm
Temperature: 0 - 50°C
Waterproof: IP57

HOT SALE

800.-

Spray (Autoclavable) 500mL
Item no. 2-7625-01
Disinfection ethanol and reverse soap can be used.

VIOLAMO

Centrifuge Tube

Centrifugal strength 20,000G !!

- Maximum centrifugal strength of 20,000 G
- DNase/RNase/DNA/progen-free
- Temperature: -80 to +121 °C
- Size: 15mL / φ17×118mm, 50mL / φ30×114mm
- Electron beam sterilized

Volume	Price/ case	Price/ 2 cases	Price/ 3 Case
15 ml Item no. 1-3500-21 50 pieces/bag x 8 bags (400 pcs)	3,800.-	7,400.-	10,000.-
50 ml Item no. 1-3500-22 20 pieces/bag x 10 bags (200 pcs)	3,200.-	6,200.-	9,000.-

หมายเหตุ: 1. สินค้ายังไม่รวมภาษี 7% 2. สินค้าไฟ 100V แคมพร้าหม้อแปลงไฟ





2026 PROMOTION

1 JANUARY - 31 DECEMBER

LAB EQUIPMENTS & CONSUMABLE



pH meter model AB23pH W/Electrode

pH



19,500.-

pH range: 0.00 to 14.00 pH
Resolution: 0.01 pH
Accuracy: +/- 0.1 pH
ORP: +/- 1999 mV
Temp: 0.0 to 100.0 °C
Display: 5" segment LCD
pH Buffer Mini Kits (4 x 50mL)

pH meter model AB41pH W/Electrode and Stirrer Compact

pH



47,000.-

pH range: -2.000 to 20.000 pH
Resolution: 0.1/0.01/0.001 pH
Accuracy: +/- 0.002 pH
ORP: +/- 2000.0 mV
Temp: -10.0 to 125.0 °C
Display: 6.5" LCD Backlight
pH Buffer Mini Kits (4 x 50mL)

EC meter model AB23EC W/Electrode

EC



29,900.-

Conductivity / TDS / Salinity
EC range: 00.0 µS/cm – 199.9 mS/cm
Resolution: 0.01 µS/cm
Accuracy: ± 1 % Reading ± 3 LSD
Display: LCD with backlight
Conductivity Standard Kits (1413 us)

Guardian™ 2000 Hotplates & Stirrers Model e-G21HS07C



11,500.-

Speed: 200 to 1600 rpm
Control: Analog
Capacity: 15 L
Temp.: 70 to 500 °C

Vortex Mixer model VXMNAL



9,990.-

Speed: 300 to 2800 rpm
Touch mode operation
Continuous mode

1 Block Dry Block Heaters Model HB1DG



Get Free 1 block

35,500.-

Control: Digital
Temp. Range: Ambient + 5°C – 120°C
Timer: 1 second - 9999 minutes
Stability: ± 0.2°C

2 Block Dry Block Heaters Model HB2DG



Get Free 2 blocks

48,500.-

Control: Digital
Temp. Range: Ambient + 5°C – 120°C
Timer: 1 second - 9999 minutes
Stability: ± 0.2°C

Microcentrifuge Tube Blocks

Select Block

Brand/Sample Type	Number of Wells	Wells Diameter	Well Depth	Item Number
0.5 mL Tube	30	7.9 mm	27.6 mm	30400157
1.5 mL Tube	20	11.1 mm	39.1 mm	30400159
1.5 mL Eppendorf Tube	20	11.5 mm	36.9 mm	30400162
2 mL Eppendorf Tube	20	11.1 mm	38.1 mm	30400191
2 mL Corning Tube	20	10.9 mm	38.1 mm	30400192



หมายเหตุ: สินค้ายังไม่รวมภาษี 7%





2026

PROMOTION

1 JANUARY - 31 DECEMBER

LAB EQUIPMENTS
& CONSUMABLE

Cole-Parmer®

**130-Watt Ultrasonic Processor with
Timer and Pulsar Item no.04714-51**



240,000.-

Sample Volume: 150 uL to 150 ml
Timer controls: 1s to 10hr.
Included: Probe with tip

**750-Watt Ultrasonic Homogenizer
with Temp Controller
Item no.04711-45**



แถมฟรี!
พร้อม Temperature probe
& Stand & Clamp
& Lab jack

378,000.-

Sample Volume: 250 uL to 1000 ml
Timer controls: 1s to 10hr.
Included: Probe with tip

**Vacuum Aspiration System,
Barbed Tube Fittings Cap, 2 L
Item no.13050-51**



34,000.-

Adjustable vacuum level, 1 mL/sec to
15 mL/sec flow rate
Max Temperature (° C): 40

**Digital Colony Counter with Magnifier
Item no. 14212-04+14212-61**



แถมฟรี!
แว่นขยาย
3.0X

49,990.-

Accept with 50- to 90-mm Petri dishes
bright-white, LED light source
Magnifier 3.0xDIGITAL

**VapLock™ Solvent Waste Kit with
Square Carboy
Item no. 12017-94**



16,700.-

10-Port, Polypropylene and HDPE, 63 mm
Closure,Capacity 10 L, include all the required
components for safe waste containment, Use the
waste kits to minimize leaks, spills, and vapor escape
from solvent waste containers

**VapLock™ Exhaust Filter
with Indicator, 75 g activated carbon
Item no.12019-80**



4,500.-

**HEATING MANTLE, 500 ML
MODEL HM-200-500
Item no. 36002-14**



22,900.-

500 mL flask capacity
Max temperature (° C): 450
Chemical resistant PP housing
is "safe to the touch"

**HEATING MANTLE, 1000 ML
MODEL HM-200-1000
Item no. 36002-18**



24,000.-

1000 mL flask capacity
Max temperature (° C): 450
Chemical resistant PP housing
is "safe to the touch"

**GAST DOA-P504-BN VACUUM PUMP
Item no.07061-42**



แถมฟรี!
Air vent filter
+ สายยาง 1 m

35,500.-

Max Flow Rate (L/min): 28
Max Vacuum (in Hg): 24
Motor hp: 1/8

**Pipette Controller OmegaZen
(Only Stock)
Item no. 25300-96**



13,500.-

Accepts glass or plastic serological
pipette 1 to 100 mL
Silicone adapter autoclavable
LCD Display

**pH buffer pH 4, 7, 10
(500 ml/bottles)**



950.-/ bottle

**Electrode Storage Solution
Item no. 00653-04
Capacity: 500 mL**



1,790.-



2026


PROMOTION

1 JANUARY - 31 DECEMBER

LAB EQUIPMENTS
& CONSUMABLE

Cole-Parmer®

**Cole-Parmer® Irreversible High-Temperature
6-Point Horizontal Labels / Indicators
(29 to 290 °C)**



- Used as a quick inspection tool that the product is in usable condition during the shipping and receiving process.
- Save time and money due to costly recalls.
- Self-adhesive design.
- EXPRESS SHIPPING BY AIR IS REQUIRED TO PREVENT ACTIVATION WHILE IN TRANSIT



No.	Item no.	Temperature points (°C)	Max. Temp (°C)
1	90308-91	29, 34, 37, 40, 42, 44	44
2	90308-92	44, 46, 49, 54, 60, 62	62
3	90308-93	65, 71, 77, 82, 88, 93	93
4	90308-94	99, 104, 110, 116, 121, 127	127
5	90308-95	132, 138, 143, 149, 154, 160	160
6	90308-96	166, 171, 177, 182, 188, 193	193
7	90308-97	199, 204, 210, 216, 224, 232	232
8	90308-98	241, 249, 254, 260, 280, 290	290

PRICE
1,700.-/pk

*Item no./10 PCS/PK



**3,900.-
/Unit**



Large and clear 4 digit display Locked -Position



**VWR® PHP Pipettes are fully autoclavable,
High performance pipettes with low plunger force.**

- Durable tip ejector made using special corrosion resistant alloy, protecting against chemical & physical corrosion
- Effortless tip ejection with unique shock absorbing mechanism
- Fully Autoclavable
- Ergonomic design with proper weight distribution
- Calibrated as per ISO 8655-6 in ISO 17025 accredited lab
- Color coding for easy identification
- Easy cleaning & in-house calibration
- Highly durable PVDF tip cone

GET FREE

**15,900.-
/Set**



Pipette Stand

Item no.	Description
VWRA613-6688	Single Channel Variable Volume 0.1-2.5 µl
VWRA613-6689	Single Channel Variable Volume 0.5-10 µl
VWRA613-6690	Single Channel Variable Volume 2-20 µl
VWRA613-6692	Single Channel Variable Volume 10-100 µl
VWRA613-6693	Single Channel Variable Volume 20-200 µl
VWRA613-6694	Single Channel Variable Volume 100-1000 µl
VWRA613-6695	Single Channel Variable Volume 1000-5000 µl
VWRA613-6696	Single Channel Variable Volume 1000-10000 µl

*คลิกขนาดได้

หมายเหตุ: สินค้ายังไม่รวมภาษี 7%

2026

PROMOTION

1 JANUARY - 31 DECEMBER



LINE Official Account

LAB EQUIPMENTS & CONSUMABLE



Photometer model DEN-600



- Estimated total number of cells (OD600)
- Turbidity measurement method (McFarland(McF))
- Protein concentration measurement (Bradford protein assay method)
- The cell holder accommodates standard 10 mm path cuvettes, round bottom, conical vials or falcon tubes.

57,000.-

Heated Circulators model T100 with T-Clamp



- Cooling/heating range: 0°C* to 100°C (*with Cooling unit)
- Simple, intuitive user interface: dial and two function buttons
- Clear, bright 4 digit LED display
- Visual alarm and low-liquid detection

78,000.-

Heated Circulators with Bath model T100-ST18



- 18l stainless steel tank
- Cooling/heating range: 0°C* to 100°C (*with Cooling unit)
- Simple, intuitive user interface: dial and two function buttons
- Clear, bright 4 digit LED display
- Visual alarm and low-liquid detection

115,700.-

Heated Circulators model TX150 with T-Clamp



- Cooling/heating range -15°C* to 150°C (*with Cooling unit)
- Icon driven home screen via dial and two function buttons
- Large, bright full colour display
- USB/RS232 interface for programming or data logging
- External probe socket
- 5 point user calibration

126,700.-

Heated Circulators with Bath model TX150-ST18



- 18l stainless steel tank
- Cooling/heating range -15°C* to 150°C (*with Cooling unit)
- Icon driven home screen via dial and two function buttons
- Large, bright full colour display
- USB/RS232 interface for programming or data logging
- External probe socket
- 5 point user calibration

167,000.-



LabT / UD15



LabT / MiniPump

26,000.-

แถมฟรี! สายยาง 2m, เลือก ID ขนาดใดขนาดหนึ่ง

- Control method: Digital knob and Mechanical keypad
- Keypad lifetime: 300,000 times
- Display: OLED display(0.96")
- Dimension: 180x115x125mm(LxWxH), 800g
- Condition temperature: 0-40°C
- IP rate: IP31

Pump Head	Speed (rpm)	Tubing Size (ID x Wall Thickness, mm)	Flow Rate (ml/min)
UD15	0.1 - 350	16#, 17#, 25#	0.08 - 930
MiniPump01	0.1 - 300	13#, 14#, 16#, 25#, 1x1, 2x1, 2.5x1, 3x1	0.0024 - 190
MiniPump02	0.1 - 300	1x1, 2x1, 2.5x1, 3x1	0.005 - 108.39



BenchMate C8-M Micro Centrifuges



- Rotor Capacity: 8x1.5/2ml Tube rotor & 2 Place PCR Strip Rotor
- Speed: 6000 RPM, 2000xg (MAX RCF)

7,900.-

BenchMate C6VX Clinical Centrifuge



- Fixed Angle Rotor Capacity: 8 x 10/15ml Adaptors for smaller 5ml, 7ml, 9/10 ml & 1.5/2 ml tubes
- Speed Setting: 500 to 6500 RPM
- Run time: 1 min to 30 mins & infinite mode

44,900.-

BenchMate C12-V High Speed Micro centrifuge



- Rotor 12-place 1.5 mL /2.0 mL microcentrifuge tube
- Maximum rpm/RCF: 15,000 rpm/15,596 x g
- Timer Setting: Countdown Timer from 1 Minute to 999 Minutes

46,900.-

Accuhelp Pipette Controller



- Volumetric Range: Plastic or Glass pipette from 0.1 mL up to 100 mL
- Battery usage Duration: Eight hours continuous use
- Choose one color: Blue, Green, Orange and Red

9,500.-/Unit

หมายเหตุ: สินค้ายังไม่รวมภาษี 7%

2026

PROMOTION

1 JANUARY - 31 DECEMBER



LINE Official Account

LAB EQUIPMENTS
& CONSUMABLE

accuhance
Enhance Your Lab Experience



- Made of ultra-pure polystyrene material that is in line with USP VI standards.
- High visibility and accurate graduation ensures easy reading of pipetting volume.
- Colour-coded for ease of identification.
- Filter barrier to prevent liquid overflow.
- Designed with two-way and negative graduation line.
- Individually packed to ensure safety and sterility.

Part No.	Volume	Packaging Size	Colour Code	Price
AC-SERO-005	5ml	50 pcs x 4 pk/case (200pcs/case)		1,550.-
AC-SERO-010	10ml	50 pcs x 4 pk/case (200pcs/case)		1,800.-
AC-SERO-025	25ml	25 pcs x 8 pk/case (200pcs/case)		2,300.-
AC-SERO-050	50ml	25 pcs x 4 pk/case (100pcs/case)		2,600.-

GVS

SYRINGE FILTER



Non-Sterile



Sterile



HOT PRICE

Material	Pore Size	Diameter	Non-Sterile (100pcs/box)	Sterile (50pcs/box)
Nylon	0.2 um	13 mm	1,700.-	-
		25 mm	1,820.-	2,100.-
	0.45 um	13 mm	1,700.-	-
		25 mm	1,820.-	2,000.-
PTFE	0.2 um	13 mm	1,890.-	-
		25 mm	2,290.-	-
	0.45 um	13 mm	1,890.-	-
		25 mm	2,250.-	-
PES	0.2 um	13 mm	1,710.-	1,800.-
		25 mm	1,970.-	2,000.-
	0.45 um	13 mm	1,710.-	1,800.-
		25 mm	1,970.-	2,000.-
CA	5.0 um	25 mm	2,250.-	-

VACUUM PUMP



19,000.-

Vacuum Pump
Model V400
Max Pressure: 100 mbar
Max flow rate: 34 L/min

AIR VENT FILTER



380.-

Air vent filters (Non-sterile)
63 mm. (dia.)
Membrane: PTFE
Pore Size: 0.2 µm

MICROCENTRIFUGE TUBE



1.5 ml



2.0 ml



5.0 ml

Volume	Price Pack	
	1-10 PK	≥ 11 pk
1.5 ml (500 pcs/pk)	380.-	300.-
2.0 ml (500 pcs/pk)	550.-	490.-
5.0 ml (200 pcs/pk)	450.-	430.-

Hycon
Plastics

PS Sterilized Petri
Dish Plastic
90×15 mm



QTY. Case	Price 500 pairs/case
1-10	1,790.-
≥ 11	1,700.-

หมายเหตุ: สินค้ายังไม่รวมภาษี 7%

ALC



Pipette Tips (China)
10 uL, 1000 ea/pk
200 uL, 1000 ea/pk
1000 uL, 500 ea/pk

BUY 3
GET 1
FREE

250.-/pk