



SampleScan Plus

Single Tube 1D & 2D Barcode Scanner




Features

- Scans and decodes in less than one second
- Improved scan alert sound and visible LED to confirm scan complete. Now louder!
- Compatible with 2D data matrix barcodes
- Compatible with 1D barcodes such as code 128, code 39, and more
- Plug-and-play barcode reader for 1D and 2D barcodes
- Keyboard wedge functionality

Labware Compatibility

Compatible with all 2D barcoded tubes from BioMicroLab, Biosero, Corning, FluidX, Greiner, Micronic, Thermo Scientific (Matrix, ABgene, Nunc), Ziath; 4mL glass vials, 2mL Cryovials and more.



Outputs scanned barcodes directly to any open application

Specifications:	Electrical	System Requirements	Dimensions	Weight
	One USB connection No power supply needed	PC and Mac Compatible	4" D X 4" W X 5.5" H 10 cm x 10 cm x 14 cm	2.5 lbs 1.13 kg



BioMicroLab provides automation solutions for tube and vial handling, liquid handling, capping/decapping, labeling, volume detection, and 2D barcode readers for pharmaceutical, biotech, agricultural and research laboratories worldwide. We have a no-compromise commitment to equipping our customers with innovative and world-class sample management tools.