

Cole-Parmer® Carboys

Get the right material for your application

- Easy-open polypropylene (PP) closure provides leakproof top connections
- Simplified transporting—full-grip handles on top and bottom make it easy to lift and carry
- Chemical-resistant PP is BPA-free, heat-stable and autoclavable
- Chemical-resistant HDPE is ideal for buffers, reagents, solids and powders
- Use amber HDPE for light- or UV-sensitive products
- Transparent PET is virtually shatterproof and withstands cold temperatures
- Some models feature an autoclavable PP spigot
- Eliminate contamination with C-Flex® or Sani-Tech® high-purity closures; call for details

NEW



Liters (gallons)	Cap size	Spigot	PP	HDPE	Amber HDPE	PET
			Catalog number	Catalog number	Catalog number	Catalog number
2.5 (0.7)	83 mm	No	JB-06059-49	JB-06059-37	JB-06058-73	JB-06059-64
		Yes	JB-06059-50	JB-06059-38	JB-06058-74	JB-06059-65
5 (1.3)	83 mm	No	JB-06059-51	JB-06059-39	JB-06058-75	—
		Yes	JB-06059-52	JB-06059-40	JB-06058-76	JB-06059-66
10 (2.6)	83 mm	No	JB-06059-53	JB-06059-41	JB-06058-77	—
		Yes	JB-06059-54	JB-06059-42	—	—
20 (5.3)	83 mm	No	JB-06059-55	JB-06059-43	—	—
		Yes	JB-06059-58	JB-06059-44	—	—

JB-06059-69 Stainless steel dolly to support 40- and 75-L carboys, max 90 kg



VapLock™ Solvent Collection Systems with Quick Disconnects and FM-Approved Container

Complete solution to safely collect solvent waste from laboratory instruments

- FM-approved container offers added safety
- Collect solvent waste from up to 10 effluent streams
- Minimize exposure to harmful chemical vapors

NEW

Systems come with an FM-approved, OSHA, NFPA, and EPA compliant 20-L polyethylene (PE) safety container with quick-disconnect connections plus multiple safety features. A pressure relief lid automatically vents between 3 and 5 psi (0.2 and 0.35 bar) reducing any risk of fire, and a stainless steel arrester prevents flashback ignition. An integral valve prevents vapor release and leaks. A large, self-closing opening and an easy-carry handle with fixed-position latch keeps lid open for convenient pouring. Select from PP or PTFE manifold with quick disconnect fitting to connects to 10 waste lines.

Material	Catalog number
PP	JB-12020-69
PTFE	JB-12020-70

